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Referenzen

- Coronaviridae Study Group of the International Committee on Taxonomy of Viruses. The species severe acute respiratory syndrome-related coronavirus: classifying 2019-nCoV and naming it SARS-CoV-2 [J]. Nature Microbiology, 5, 536-544 (2020).
- Perlman, S. Netland, J. Coronaviruses post-SARS: update on replication and pathogenesis. Nature Reviews Microbiology 7, 439-450, doi: 10.1038/nrmicro2147 (2009).
- Lauer SA, Grantz KH, Bi Q, et al. The Incubation Period of Coronavirus Disease 2019 (COVID-19) From Publicly Reported Confirmed Cases: Estimation and Application. Ann Intern Med. 2020; 172 (9): 577-582. doi: 10.7326/M20-0504.
- McCarron, MM, et al. Detection of Phencyclidine Usage by Radioimmunoassay of Saliva. J Anal Toxic.1984 Sep-Oct.; 8 (5), pp 197-201.
- B. Korber et al. Tracking Changes in SARS-CoV-2 Spike: Evidence that D614G Increases Infectivity of the COVID-19 Virus. Cell, vol. 182, no. 4, pp. 812-827.e19, Aug. 2020, doi: 10.1016/j.cell.2020.06.043.
- Investigation of novel SARS-CoV-2 variant: Variant of Concern 202012/01, GOV.UK. https://www.gov.uk/government/publications/investigation-of-novel-sars-cov-2-variant-variant-of-concern-20201201 (accessed Jan. 19, 2021).
- New SARS-CoV-2 variant, GOV.UK. https://www.gov.uk/government/collections/new-sarscov-2-variant (accessed Jan. 11, 2021).
- D. C. Dinesh et al. Structural basis of RNA recognition by the SARS-CoV-2 nucleocapsid phosphoprotein. PLoS Pathog, vol. 16, no. 12, p.e1009100, Dec. 2020, doi:10.1371/journal.ppat.1009100.
- J. Mariën et al. Evaluating SARS-CoV-2 spike and nucleocapsid proteins as targets for antibody detection in severe and mild COVID-19 cases using a Lumindex bead-based assay. Journal of Virological Methods, vol. 288, p. 114025, Feb. 2021, doi: 10.1016/j.jviromet.2020.114025.

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Instructions for use

Intended use and Principle

IME® SARS-CoV-2 Ag Saliva Rapid Test is for the rapid, qualitative detection of the nucleocapsid protein antigen from SARS-CoV-2 in human saliva^[1]. The test is for *in vitro* diagnostic use only. For professional use only. It is intended for clinical laboratories and healthcare professional use only for point-of-care testing^[2]. It provides only an initial screening test result. More specific alternative diagnosis methods (molecular diagnostic and/or CT) should be performed in order to obtain the confirmation of SARS-CoV-2 infection. Not for at-home testing.

IME® SARS-CoV-2 Ag Saliva Rapid Test is based on immunoassay technology. Each test device has one line of anti-SARS-CoV-2 monoclonal antibody on the detection line (T line) and one line of anti-mouse IgG polyclonal antibody on the quality control line (C line). When extracted specimen is added to the specimen well, it will react with the labeled antibody to form a complex, the mixture then migrates through the membrane by capillary action and interacts with the coated anti-SARS-CoV-2 monoclonal antibody on the detection line. If the specimen contains SARS-CoV-2 antigen, the detection line will appear purplish-red indicating the SARS-CoV-2 antigen is positive. Otherwise, the test result will be negative. The test device also contains a quality control line C which should appear purplish red for all valid tests. If the quality control line C does not appear, the test result will be invalid even if the detection line appears.

Package Content

Each test kit contains test devices, sealed pouches (prefilled with extraction solution), extraction tubes, extraction tube tips, collection cups, pipettes (fixed volume), tube stand and package insert. *Materials required but not provided: timer.*

Storage and Handling

- Store the test kit in a cool, dry place between 2 °C – 30 °C. Keep away from light.
- Exposure to temperature and/or humidity outside the specified conditions may cause inaccurate results.
- Do not freeze. Use the test kit at temperatures between 15 °C – 30°C.
- Use the test kit between 10 % - 90 % humidity.
- Do not use the test kit beyond the expiration date (printed on the foil pouch and box). **Note:** All expiration dates are printed in Year-Month-Day format. 2022-06-18 indicates June 18, 2022.


Warning, Precautions and Limitations

- Results from SARS-CoV-2 antigen testing should not be used as the sole basis to diagnose or exclude SARS-CoV-2 infection or to inform infection status^[2].
- Negative results do not rule out SARS-CoV-2 infection, particularly in those who have been in contact with the virus. Follow-up testing with a molecular diagnostic and/or CT should be considered to rule out infection in these individuals.
- Positive results may be due to present infection with SARS-coronavirus strains, see **Cross-Reactivity** for details. Follow-up testing with a molecular diagnostic and/or CT should be considered to confirm the testing result.
- False positive results may be due to sticky specimen, insufficient specimen volume or bubbles during applying.
- Please take a new pipette to draw the specimen if the pipette is damaged or cannot use.
- Do not use unverified UTM, which may lead to false positive or false negative results.
- For *in vitro* diagnostic use only. Not for at-home testing.
- Further molecular diagnostic and/or CT is recommended to identify the actual physical situation.
- Do not open the foil pouch of the test device exposing it to the ambient environment until the test device is ready for immediate use.
- Do not use any damaged test device or material.
- Do not reuse the test device.
- Handle extraction solution with caution, do not contact with eyes or skin. If spilled on eyes or skin, wash thoroughly with water.
- Do not use test kit beyond the expiration date.
- Wash hands thoroughly after handling.
- Specific training or guidance is recommended if operators are not experienced with specimen collection and handling procedures.
- Only use saliva as specimen. Follow the package insert to obtain accurate results.
- Wear protective gears such as laboratory coats, disposable gloves, and eye protection when specimens are collected and evaluated.
- All parts of kit are considered biohazardous and can potentially transmit infectious diseases from pathogens, even after you have performed cleaning and disinfection. Follow proper precautions and all local regulations when disposing of the used test kits.
- Since the outbreak of the pandemic, the SARS-CoV-2 variant with D614G mutations in the spike protein has replaced the original form in the most regions all over the world^[3]. In December 2020, a novel strain of the virus, named 'VUI-202012/01', was identified in the England with a set of 17 mutations^[4]. Another mutant strain 501Y.V2 of SARS-CoV-2, originally detected in South Africa, shares the same key mutation N501Y. The N501Y mutation locates the receptor-binding domain (RBD) of the spike protein that the virus uses to bind to the human ACE2 receptor, which might be associated with increased transmissibility^[7]. The nucleocapsid phosphoprotein (N Protein), linking the viral envelope to the viral RNA, plays a central role in the packaging signal RNA recognition and subsequent RNA encapsidation^[8]. Based on its vital role in transcription and replication of the virus, the N protein is suggested to be more sensitive for the early detection of infections^[9]. The present SARS-CoV-2 Ag rapid test employs the interaction with antigen sites in N protein. Till now, there is no clear evidence indicating that mutations found in Spike protein can affect the performance of N protein based antigen tests.

Specimen Collection^[4]

Please instruct the person who need to be tested not to place anything in the mouth including food, drink, gum, tobacco, etc. for at least 10 minutes prior to collection.

1. Remove a collection cup from the box.
2. Spit the saliva into the collection cup (1).
3. Fill the pipette from the collection cup with saliva:
 - 3.1. Firmly squeeze the top bulb.
 - 3.2. Still squeezing, place the pipette tip into the saliva.
 - 3.3. With the pipette tip still in the saliva, slowly release pressure on bulb to fill the pipette until the saliva fills the overflow bulb (2).

Note:

- Stored at room temperature or cooled down (2 °C - 8 °C), the samples can be used for 24 hours without influencing the test results.
- Do not freeze the specimen.
- Please recollect the specimen if the volume is less than 500 µL. Make sure the specimen volume is more than 500 µL.
- Do not use any damaged material.

Specimen Handling

Freshly collected specimens should be tested as soon as possible (not more than 24 hours after the specimen is collected into the collection cup). It is essential that correct specimen collection and preparation methods are followed.


Test Procedure

Allow the Test Devices and Extraction Solution to equilibrate to 15 °C – 30 °C prior to testing.

1. Hold the sealed pouch vertically and let all extraction solution flow into the bulb. Break the tip (3) and squeeze the bulb to dispense all extraction solution into the extraction tube (4).
2. Collect specimen refer to **Specimen Collection**.
3. Firmly squeeze the top bulb to empty the contents of pipette into the extraction tube with extraction solution (5). Extra liquid over in the overflow bulb should be left behind. Shake the tube 10 times (6). Please shake the tube slowly and gentle to avoid the bubbles.

Note: Do not use the extraction solution if it leaks. Use a new sealed pouch (prefilled with extraction solution) and follow steps 1 to 3 to extract the specimen. Please mix specimen well with the extraction solution if the specimen is sticky.
4. Put on the tube tip (7). Please proceed with the extracted specimen within the next 60 minutes.

5. Take out a test device from sealed foil pouch and put it on a clean and level surface.
6. Apply 3 drops of the extracted specimen into the specimen well by pressing in the middle of the extraction tube smoothly (8). Please avoid bubbles during application.
7. Read the test result at 15 minutes (9). Don't read the result after 20 minutes.


Note:

- Do not interchange or mix extraction solution from different lots.
- Please follow local regulations to handle the used materials.

Interpretation of Test Results

1. **Positive Result:** Both the quality control line C and the detection line T appear (10).
2. **Negative Result:** Only the quality control line C appears, with no other line appearing on the detection line (11).
3. **Invalid Result:** Quality control line C fails to appear indicating the test is invalid, no matter if the detection line T appears or not (12). In this case collect a new specimen and perform another test with a new test device.

Quality Control

Internal procedural controls are included in the test. A colored line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

Performance

I. Limit of Detection

The LOD for the IME® SARS-CoV-2 Ag Saliva Rapid Test is established using dilutions of an inactivated virus culture. The starting material is supplied at a concentration of 1.51 x 10⁶ TCID₅₀/mL. Studies are designed to estimate the LOD of the assay using saliva specimens, the starting material is spiked into a volume of pooled human saliva matrix obtained from healthy donors and confirmed negative for SARS-CoV-2 to obtain a series of different concentrations.

SARS-CoV-2 Titer	1.51 x 10 ⁶ TCID ₅₀ /mL							
Dilution	1/10	1/100	1/1000	1/2500	1/5000	1/10000	1/20000	1/40000
Concentration in Dilution tested (TCID ₅₀ /mL)	1.51 x 10 ⁵	1.51 x 10 ⁴	1.51 x 10 ³	6.04 x 10 ²	3.02 x 10 ²	1.51 x 10 ²	75.5	37.8
Detection rates of 5 replicates	100 % (5/5)	100 % (5/5)	100 % (5/5)	100 % (5/5)	100 % (5/5)	100 % (5/5)	100 % (5/5)	60 % (3/5)
Detection rates of 20 replicates near cut off	NA	NA	NA	NA	100 % (20/20)	100 % (20/20)	95 % (19/20)	75 % (15/20)
Lowest Concentration with Uniform Positivity per Analyte	75.5 TCID ₅₀ /mL							
Limit of detection (LoD) per inactivated Virus Culture	75.5 TCID ₅₀ /mL							

II. Clinical Sensitivity and Clinical Specificity

A total of 603 specimens are tested using the IME® SARS-CoV-2 Ag Saliva Rapid Test. These specimens are obtained by saliva from symptomatic patients. The performance of the IME® SARS-CoV-2 Ag Saliva Rapid Test is compared to a commercialized molecular assay.

IME® SARS-CoV-2 Ag Saliva Rapid Test	PCR		
	Positive	Negative	Total
Positive	169	0	169
Negative	3	431	434
Total	172	431	603
Sensitivity	98.26 % (169/172 95% CI 95 % - 99.41 %)		
Specificity	> 99.99 % (431/431 95 % CI 99.12 % - 100 %)		
Accuracy	99.50 % (600/603 95 % CI 98.55 % - 99.83 %)		

▲ Table Summary of sensitivity/specificity of the Ag Saliva Rapid Test compared to PCR.

The IME® SARS-CoV-2 Ag Saliva Rapid Test shows a:

- clinical sensitivity of 98.26 %.
- clinical specificity of > 99.99 %.
- clinical accuracy of 99.50 %.

Note:

Results from SARS-CoV-2 antigen testing should not be used as the sole basis to diagnose or exclude SARS-CoV-2 infection or to inform infection status. Follow-up testing with a molecular diagnostic and/or CT should be considered to confirm the testing result.

Cross-Reactivity

I. Cross-Reactivity

There is no cross-reaction with potential cross-reactive substances except SARS-Coronavirus.

Cross-reaction with SARS-coronavirus:

Virus	Strain	Concentration
SARS-coronavirus	urbani	1 x 10 ⁶ PFU/mL

No cross-reaction with potential cross-reactive substances:

Virus/Bacteria/Parasite	Strain	Concentration Range
Influenza A	H1N1	1 x 10 ⁴ ~ 1 x 10 ⁶ TCID ₅₀ /mL
	H3N2	
	H5N1	
	H7N9	
	N/A	
Influenza B	Type 1	1 x 10 ⁴ ~ 1 x 10 ⁶ TCID ₅₀ /mL
	Type 2	
	Type 3	
	Type 5	
Adenovirus	Type 7	1 x 10 ⁴ ~ 1 x 10 ⁶ TCID ₅₀ /mL
	Type 55	
	Type B	
Respiratory syncytial virus	Type A	1 x 10 ⁴ ~ 1 x 10 ⁶ TCID ₅₀ /mL
	Type B	
Coronavirus	229E	1 x 10 ³ ng/mL
	OC43	
	NL63	
	HKU1	

Virus/Bacteria/Parasite	Strain	Concentration Range
MERS-Coronavirus	Florida/USA-2_Saudi Arabia.2014	1 x 10 ⁴ ~ 1 x 10 ⁶ TCID ₅₀ /mL
	Type 1	
	Type 2	
	Type 3	
Parainfluenza virus	Type 4	1 x 10 ⁴ ~ 1 x 10 ⁶ TCID ₅₀ /mL
	Type 4	
Rhinovirus A16	N/A	
Legionella pneumophila	Bloomington-2 82A3105	
Mycobacterium tuberculosis	K	1 x 10 ⁵ cells/mL
	Erdman	
	HN878	
	CDC1551 H37Rv	
Streptococcus pneumonia	475298 [Maryland(D1)6B-17]	1 x 10 ⁵ cells/mL
	178 [Poland23F-16]	
	262 [CIP 104340]	
Streptococcus pyrogens	Slovakia 14-10 [29055]	1 x 10 ⁵ cells/mL
	Typing stain T1	
Mycoplasma pneumoniae	Mutant22	1 x 10 ⁵ cells/mL
	FH strain of Eaton Agent M129-B7	

II. Interference Substances Studies

There is no interference for potential interfering substances listed below:

Potential Interfering Substance	Concentration	
Anti-viral drugs	Zanamivir (Influenza)	5 mg/mL
	Oseltamivir (Influenza)	10 mg/mL
	Artemether-lumefantrine (Malaria)	50 µmol/L
	Doxycycline hyclate (Malaria)	70 µmol/L
	Quinine (Malaria)	150 µmol/L
	Lamivudine (Retroviral medication)	1 mg/mL
	Ribavirin (HCV)	1 mg/mL
Respiratory Specimens	Daclatasvir (HCV)	1 mg/mL
	Mucin: bovine submaxillary gland, type I-S	100 µg/mL
	Blood (human), EDTA anticoagulated	5 % (v/v)
Nasal sprays or drops	Biotin	100 µg/mL
	Neo-Synephrine (Phenylephrine)	10 % (v/v)
	Afrin Nasal Spray (Oxymetazoline)	10 % (v/v)
Homeopathic allergy relief medicine	Saline Nasal Spray	10 % (v/v)
	Homeopathic Zicam Allergy Relief Nasal Gel	5 % (v/v)
	Sodium Cromoglycate	20 mg/mL
Anti-inflammatory medication	Olopatadine Hydrochloride	10 mg/mL
	Acetaminophen	199 µg/mL
	Acetylsalicylic acid	3.62 mmol/L
Antibiotic	Ibuprofen	2.425 mmol/L
	Mupirocin	10 mg/mL
	Tobramycin	5 µg/mL
	Erythromycin	81.6 µmol/L
	Ciprofloxacin	30.2 µmol/L

III. High-Dose Hook Effect

Cultured SARS-CoV-2 virus is spiked into specimen. No hook-effect was observed at 1.51 x 10⁶ TCID₅₀/mL of cultured SARS-COV-2 virus.

References

- [1] Coronaviridae Study Group of the International Committee on Taxonomy of Viruses. The species severe acute respiratory syndrome-related coronavirus: classifying 2019-nCoV and naming it SARS-CoV-2 [J]. Nature Microbiology, 5, 536-544 (2020).
- [2] Perlman, S. Netland, J. Coronaviruses post-SARS: update on replication and pathogenesis. Nature Reviews Microbiology 7, 439-450, doi: 10.1038/nrmicro2147 (2009).
- [3] Lauer SA, Grantz KH, Bi Q, et al. The Incubation Period of Coronavirus Disease 2019 (COVID-19) From Publicly Reported Confirmed Cases: Estimation and Application. Ann Intern Med. 2020; 172 (9): 577-582. doi: 10.7326/M20-0504.
- [4] McCarron, MM, et al. Detection of Phencyclidine Usage by Radioimmunoassay of Saliva. J Anal Tox.1984 Sep-Oct.; 8 (5), pp 197-201.
- [5] B. Korber et al. Tracking Changes in SARS-CoV-2 Spike: Evidence that D614G Increases Infectivity of the COVID-19 Virus. Cell, vol. 182, no. 4, pp. 812-827.e19, Aug. 2020, doi: 10.1016/j.cell.2020.06.043.
- [6] Investigation of novel SARS-CoV-2 variant: Variant of Concern 202012/01, GOV.UK. https://www.gov.uk/government/publications/investigation-of-novel-sars-cov-2-variant-variant-of-concern-20201201 (accessed Jan. 19, 2021).
- [7] New SARS-CoV-2 variant, GOV.UK. https://www.gov.uk/government/collections/new-sarscov-2-variant (accessed Jan. 11, 2021).
- [8] D. C. Dinesh et al. Structural basis of RNA recognition by the SARS-CoV-2 nucleocapsid phosphoprotein. PLoS Pathog, vol. 16, no. 12, p. e1009100, Dec. 2020, doi:10.1371/journal.ppat.1009100.
- [9] J. Mariën et al. Evaluating SARS-CoV-2 spike and nucleocapsid proteins as targets for antibody detection in severe and mild COVID-19 cases using a Luminex bead-based assay. Journal of Virological Methods, vol. 288, p. 114025, Feb. 2021, doi: 10.1016/j.jviromet.2020.114025.

EXPLANATION OF SYMBOLS			
	Catalogue number		Lot number
	Manufacturer		Do not reuse
	Date of manufacture		Contains sufficient for <n> tests
	Authorized representative in the European Community		Temperature limit
	Consult instructions for use		In vitro diagnostic medical device
	Use-by date		This product meets the requirements of the directive EG 98/79

REF: A712001DEEN-XX00X EAN: 0811900031670 | PZN: 17179701

	VivaChek Biotech (Hangzhou) Co., Ltd. Level 2, Block 2, 146 East Chaofeng Rd., Yuhang Economy Development Zone, Hangzhou, 311100, China Email: info@vivachek.com		Lotus NL B.V. Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands. Tel: +31644168999 Email: peter@lotusnl.com
--	--	--	---

Trademark Owner, Distribution and Importer

IME-DC GmbH
Fuhrmannstraße 11
95030 Hof | Germany
Tel: +49 9281 | 85 01 6-0 E-Mail: info@imedc.de
Fax: +49 9281 | 85 01 6-100 Web: www.imedc.de
Number: 1604025602
Effective date: 2021-01-28