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(2) Ji, W.; Wang, W.; Zhao, X.; Zan, J.; Li, X.: Homologous Recombination within the Spike Glycoprotein of the newly identified Coronavirus may boost Cross-species Transmission from Snake to Human. Journal of Medical Virology. 22. Januar 2020.

(3) Novel Coronavirus Primer and Probe Sequences, National Institute of Virology, Chinese Center for Disease Control and Prevention, 21. Januar 2020.

(4) Clinical Management of severe acute Respiratory Infection when novel Coronavirus (nCoV) Infection is suspected. Interim Guidance 2020.

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SARS-CoV-2 Nucleic Acid Testing Card FOR PROFESSIONAL USE ONLY

[Product Name]

Product Name: SARS-CoV-2 Nucleic Acid Testing Card
Trade Name: Puddle SARS-CoV-2 Card

[Intended Use]

This kit is used for *in vitro* qualitative detection of the N gene and ORF1ab gene of novel coronavirus SARS-CoV-2 in nasal swab samples from suspected pneumonia cases of SARS-CoV-2 infection and suspected cluster cases and other persons requiring diagnosis or differential diagnosis of SARS-CoV-2 infection.
The test results of this kit are only for clinical reference and shall not be used as the sole standard for clinical diagnosis. A comprehensive analysis of the patient's clinical manifestations and other laboratory tests is recommended.

[Testing Principle]

The assay is based on isothermal amplification method and enzyme digestion probe technology, and SARS-CoV-2 N gene and ORF1ab gene is selected as the target region, designed the probe and specific primers. A large number of target sequence copies were generated in the reaction system during the isothermal amplification. When the probe hybridizes to the complementary sequence, it is cleaved and fluorescence is emitted. Integrated Nucleic Acid Testing Device detects and analyzes fluorescence signal automatically, reporting negative, positive or invalid result. The assay includes internal control for monitoring of sample collection, processing, and amplification to reduce false negative results.

[Composition and Article Number]						
Article No. and specifications	RA10102 02-1 1 Test	RA10102 02-2 2 Tests	RA10102 02-5 5 Tests	RA10102 02-10 10 Tests	RA10102 02-20 20 Tests	RA10102 02-50 50 Tests
SARS-CoV-2 Reaction Card (piece)	1	2	5	10	20	50
Nucleic Acid Releasing Agent O1 (tube)	1	2	5	10	20	50
Disposable Sampling Swab (branch)	1	2	5	10	20	50
Waste Bag (piece)	1	2	5	10	20	50

- NOTE:**
- The above ingredients of different batches of kits shall not be used interchangeably.
 - The card holder is in the packing box of the Integrated Nucleic Acid Testing Device.

[Storage Conditions and Shelf-life]

- 2-8°C storage, valid for 3 months.
- The production date and expiration date are shown on the package label.

[Applicable Instruments]

Integrated Nucleic Acid Testing Device (IN4001)

[Sample Requirements]

Nasal Swab

[Testing Method]

The room temperature should be between 15 ~ 30°C. Please read all the instructions carefully before you begin.
Be careful of these operations will probably leads to an incorrect/invalid result.

STEP 1: Instruments Preparing

- Take out the device, power adapter, and card holder from Integrated Nucleic Acid Testing Device package.
- Put the Integrated Nucleic Acid Testing Device on a flat surface, connect the power supply, press the button in front of the instrument, to enter the warm-up process (the power light is flashing red). After 2 minutes, the warm-up is completed and in standby (the power light is blue).

STEP 2: Sample Collection

- Wash or disinfect your hands prior to perform the test procedure.
- Place the kit at room temperature (15 to 20°C) and take out components that required for one test.

NOTE: Make sure there are no damages of the SARS-CoV-2 Reaction Card bag, and no liquid leakage of the nucleic acid releasing agent O1. Do not use them if any leakage occurs.


NOTE: Before sample collection, maintain proper infection control and use recommended personal protective equipment (PPE), which includes an N95 or higher-level respirator (or face mask if a respirator is not available), eye protection, gloves, and a gown.

- Have the person being tested wipe away excess snot, if present.
- Take out the disposable sampling swab from the pouch and hold with handle end.
- Tip the head (or the person being tested) back (approximately 70 degrees).
- Gently insert the entire absorbent tip of the disposable sampling swab into one nostril until you meet resistance (approximately 1.5 to 2 cm).
- Rotate it 3 times against the inner wall of one nostril, so that the absorbent surface of the swab can be moistened.
- Slowly remove the disposable sampling swab from the first nostril.
- Repeat 5-6 steps with the same disposable sampling swab on the other nostril.

NOTE:

- Avoid scratches.
- Samples should be tested immediately after collection.


[Sample Processing]

STEP 3: Sample Processing

- Open the seal of the nucleic acid releasing agent O1 vial carefully to avoid spilling the liquid. Discard the seal into waste bag.
- Insert the sampled disposable sampling swab into the vial and make sure the absorbent tip is in the liquid. Then twist the disposable sampling swab tip against the bottom and sides of the nucleic acid releasing agent O1 vial 10 times and pinch the disposable Sampling Swab tip.
- Discard the disposable sampling swab into waste bag.
- Screw the cap.

NOTE: Be careful to avoid contact with eyes or skin by the nucleic acid releasing agent O1. If it happens unfortunately, wipe off the liquid immediately and rinse with plenty of water.

STEP 4: Sample Testing

Make sure the Integrated Nucleic Acid Testing Device is in standby (the power light is blue).

- Tear open the aluminum foil bag of one SARS-CoV-2 Reaction Card and take it out, place the SARS-CoV-2 Reaction Card on the card holder and unscrew the cap of the sample tube on the SARS-CoV-2 Reaction Card.

NOTE: The SARS-CoV-2 Reaction Card must proceed to subsequent operations as soon as possible after the aluminum foil bag has been torn and proceed to the next step immediately when the cap of the tube is unscrewed.

- Open the top cap of the nucleic acid releasing agent O1 vial from STEP 3, use one hand to stabilize the card holder, use the other hand to pour the nucleic acid releasing agent O1 solution on the tube inside wall of the SARS-CoV-2 Reaction Card between the two liquid injection lines by squeezing the nucleic acid releasing agent O1 vial wall.

[Sample Testing]

NOTE: There are two liquid injection lines marked on the SARS-CoV-2 Reaction Card sample tube. Add nucleic acid releasing agent O1 solution into the SARS-CoV-2 Reaction Card sample tube until the liquid level between the two lines.


- Place the nucleic acid releasing agent O1 vial in waste bag for disposal.
- Screw the cap of the SARS-CoV-2 Reaction Card sample tube tightly.
- Firmly press the protruding air bag on the sample tube cap of the SARS-CoV-2 Reaction Card to deform it and recess it into the tube.

NOTE: Wait for 10 seconds. DO NOT move on to the next step until the liquid fully fills all of the chamber that lyophilized reagent beads contained inside.

- Hold the chip and shake it up and down 10 times in about 5 seconds. The SARS-CoV-2 Reaction Card is ready to be tested.

NOTE: The SARS-CoV-2 Reaction Card must be tested immediately.

- Open the cabin door of Integrated Nucleic Acid Testing Device, and insert the SARS-CoV-2 Reaction Card into the device according to the direction indicated on the SARS-CoV-2 Reaction Card, and push it to the fixed position of the bottom card slot, close the cabin door.